

## V5-Tag

---

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of V5 epitope

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Background:**

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids, although it has also been used with a shorter 9 amino acid sequence.

*To place an order, please [Click HERE](#).*

**Catalog #:**AE019

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Isotype:**IgG

**Purity:**Affinity purification

For research use only.