www.neobiolab.com info@neobiolab.com 888.754.5670, +1 617.500.7103 United States 0800.088.5164, +44 020.8123.1558 United Kingdom

HCV NS5

Description: The E.coli derived recombinant protein contains the HCV NS5 immunodominant regions, amino acids 2061-2302. The protein is fused with GST at N-terminus.

Catalog #:HCPS-242

Purity:HCV-NS5 protein is >95% pure as determined by 10% PAGE (coomassie staining).

For research use only.

Purification Method:

HCV-NS5 protein was purified by proprietary chromatographic technique.

Specificty:

Immunoreactive with sera of HCV-infected individuals.

Formulation:

1.5M Urea, 25mM Tris-HCl pH 8, 50% glycerol and 0.2% Triton-X.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

HCV-NS5 antigen is suitable for ELISA and Western blots, excellent antigen for detection of HCV with minimal specificity problems.

Introduction:

HCV is a small 50nm, enveloped, single-stranded, positive sense RNAvirus in the family Flaviviridae. HCV has a high rate of replication with approximately one trillion particles produced each day in an infected individual. Due to lack of proofreading by the HCV RNA polymerase, the HCV has an exceptionally high mutation rate, a factor that may help it elude the host's immune response. Hepatitis C virus is classified into six genotypes(1-6) with several subtypes within each genotype. The preponderance and distribution of HCV genotypes varies globally. Genotype is clinically important in determining potential response to interferon-based therapy and the required duration of such therapy. Genotypes 1 and 4 are less responsive to interferon-based treatment than are the other genotypes (2, 3, 5 and 6).

Storage:

HCV NS5 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please Click HERE.





