

HEV Mosaic

Description: The e.coli derived HEV Mosaic protein contains 12 HEV immunodominant regions from ORF2 and ORF3 having a molecular mass of 38.5 kDa.

Catalog #:HEPS-242

Purity: HEV Mosaic protein is >95% pure as determined by 10% PAGE (coomassie staining).

For research use only.

Purification Method:

HEV Mosaic protein was purified by proprietary chromatographic technique.

Specificity:

Immunoreactive with sera HEV-infected individuals.

Formulation:

25mM Tris-HCl, 1mM EDTA, 0.5M urea and 50% glycerol.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

HEV Mosaic antigen is suitable for ELISA and Western blots, excellent antigen for detection of HEV with minimal specificity problems.

Introduction:

Hepatitis E virus (HEV), the major etiologic agent of enterically transmitted non-A, non-B hepatitis worldwide, is a spherical, non-enveloped, single stranded RNA virus that is approximately 32 to 34 nm in diameter. HEV belongs to a genus of HEV-like viruses (unassigned genus). HEV has a single-stranded polyadenylated RNA genome of approximately 8 kb. Based on its physicochemical properties it is presumed to be a calici-like virus.

Storage:

HEV Mosaic protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please [Click HERE](#).