www.neobiolab.com info@neobiolab.com 888.754.5670, +1 617.500.7103 United States 0800.088.5164, +44 020.8123.1558 United Kingdom

H3N2 Kiev

Description: Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Kiev/301/94 like /Johannesburg/33/94. The Influenza Virus was purified by Ultracentrifugation with 10-40 % sucrose gradient.

Catalog #:IHPS-012

For research use only.

Physical Appearance: Sterile Filtered colorless solution.

Purity: Greater than 90.0% as determined by Analysis by SDS-PAGE.

Formulation:

The H3N2 A/Kiev/301/94 solution contains 0.1M NaCl, 10mM Tris-Hcl, 1mM EDTA pH-8, 0.1% sodium azide (NaN3) and 0.005 % thimerosal.

Stability:

A/Kiev/301/94 although stable 4°C for 4 weeks, should be stored below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drµgs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

H3N2 is a subtype of the influenza A virus. Its name derives from the forms of the two kinds of proteinson the surface of its coat, hemagglutinin(H) and neuraminidase(N). H3N2 exchanges genes for internal proteins with other influenza subtypes. H3N2 has tended to dominate in prevalence over H1N1, H1N2, and influenza B. H3N2 strain descended from H2N2by antigenic shift, in which genes from multiple subtypes re-assorted to form a new virus. Both the H2N2and H3N2 strains contained genesfrom avian influenzaviruses.

To place an order, please Click HERE.





