www.neobiolab.com info@neobiolab.com 888.754.5670, +1 617.500.7103 United States 0800.088.5164, +44 020.8123.1558 United Kingdom

# Omp Pylori

Description: Omp Pylori recombinant antigen is produced in E. coli expressing the H. pylori outer membrane protein having the Mw of 23 kDa. Omp Pylori recombinant antigen is well recognized by specific IgG and IgM from H. pylori infected patients.

Catalog #:PRPS-376

For research use only.

Source: E. Coli

Physical Appearance: Sterile filtered liquid formulation.

## Amino Acid Sequence:

MLVTKLAPDFKAPAVLGNNEVDEHFELSKNLGKNGAILFFWPKDFTFVCPTEIIAFDKRVKDFQEK GFNVIGVSIDSEQVHFAWKNTPVEKGGIGQVTFPMVADITKSISRDYDVLFEEAIALRGAFLIDKNM KVRHAVINDLPLGRNADEMLRMVDALLHFEEHGEVCPAGWRKGDKGMKATHQGVAEYLKENSIK L.

Purity: Greater than 95% pure as determined by 12% PAGE (Coomassie staining).

### **Purification Method:**

Purified by proprietary chromatographic technique.

### Formulation:

The Omp Pylori recombinant protein is formulated in 1xPBS pH 7.4.

## Stability:

Omp Pylori although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

# Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Introduction:

Helicobacter pylori, a Gram-negative, microaerophilic bacterium that populates in the stomach, predominantly at the antrum. Helicobacter pylori causes a chronic inflammation of the mucoid lining of stomach and is highly related to the growth of duodenal and gastric ulcers and stomach cancer. Over 50% of the world's population harbor H. pylori in their upper gastrointestinal tract, infection is more prevalent in developing countries. Thus far, no ideal target H. pylori antigen has been developed for the diagnostic purpose. 23 kDa H. pylori outer membrane protein was identified with a good sensitivity and coverage in the diagnosis of H. pylori infection.

To place an order, please Click HERE.





